

Cardboard Freezer Boxes

I. APPLICATIONS

1. Storage and organization of any range of 0.5ml to 50ml tubes
2. Great for upright or chest freezers

II. FEATURES

1. Made from durable cardboard with variety of outer coatings: (Superior white , Plasti-coat, Cryo-color)
2. Chemically resistant to alcohols and mild organic solvents
3. Temperature range: stable from -196°C to 121°C
4. Suitable for liquid nitrogen storage

III. BOX DESIGN

1. Available in 2 inch and 3 inch sizes with either 81-well, 100-well, 36-well or 16-well dividers
2. Select from our Superior white coating, our Plasti-coat laminated coating, our Cryo-color or ID-color assorted colored boxes. (Pink, blue, yellow and green)
3. Holds 0.5ml, 1.5ml or 2.0ml microtubes, and 5ml, 15ml or 50ml large volume tubes
4. Removable lid includes a inventory list writing area
5. Cardboard dividers

IV. PACKAGING

1. Bulk, 100 boxes/case

V. DESCRIPTION

Item	Description	Packaging	Sterile
90-1536	36-well, superior white	100 pieces/case	No
90-5016	16-well, superior white	100 pieces/case	No
90-2281	2in, 81-well, plasti-coat	100 pieces/case	No
90-2200	2in, 100-well, plasti-coat	100 pieces/case	No
90-2381	3in, 81-well, plasti-coat	100 pieces/case	No
90-2300	3in, 100-well, plasti-coat	100 pieces/case	No
90-5281	2in, 81-well cryo-color, assorted	100 pieces/case	No
90-5200	2in, 100-well cryo-color, assorted	100 pieces/case	No

VI. STERILITY

1. N/A

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.





VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

IX. PRODUCT PARAMETER





Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
90-2381	Box	130.0	127.0	15.0	15.0	74	73.5	30.0
		H	I	J	K	L	-	-
		34.0	32.5	133.0	130.0	76.0	-	-







Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
90-2300	Box	130.0	127.0	12.5	12.5	74	73.5	30.0
		H	I	J	K	L	-	-
		34.0	32.5	133.0	130.0	76.0	-	-







Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
90-2281	Box	130.0	127.0	15.0	15.0	45.5	40.0	30.0
90-5281		H	I	J	K	L	-	-
		34.0	32.5	133.0	130.0	47.0	-	-







Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
90-2200	Box	130.0	127.0	12.5	12.5	45.5	40.0	30.0
90-5200		H	I	J	K	L	-	-
		34.0	32.5	133.0	130.0	47.0	-	-



Dimensions=mm


Cat #	Description	A	B	C	D	E	F	G
90-1536	Box	145.0	142.0	25.0	25.0	75.0	73.5	70.0
		H	I	J	K	L	-	-
		68.5	67.5	148.0	145.0	76.0	-	-



90-1536

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
90-5016	Box	145.0	142.0	35.0	35.0	75.0	73.5	70.0
		H	I	J	K	L	-	-
		68.5	67.5	148.0	145.0	76.0	-	-



90-5016